



RABGAP1L Antibody

Product Code	CSB-PA019232GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7ZAP0
Immunogen	Human RABGAP1L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	RAB GTPase activating protein 1-like;RABGAP1L;DKFZp686E1450;HHL;TBC1D18 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RABGAP1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.