



RAC3 Antibody

Product Code	CSB-PA019249GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60763
Immunogen	Human RAC3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ras-related C3 botulinum toxin substrate 3 (rho family, small GTP binding protein Rac3);
Product Type	Purified Rabbit Anti Human, Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAC3
Target Details	This protein is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases.