



RAPGEFL1 Antibody

Product Code	CSB-PA019335GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHV5
Immunogen	Human RAPGEFL1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Rap guanine nucleotide exchange factor (GEF)-like 1;RAPGEFL1;Link-GEFII;MGC134798;MGC134799 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAPGEFL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.