



RASGEF1A Antibody

Product Code	CSB-PA019357GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9B8
Immunogen	Human RASGEF1A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	RasGEF domain family, member 1A;RASGEF1A;CG4853;FLJ37817 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RASGEF1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.