



RBM14 Antibody

Product Code	CSB-PA019404GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PK6
Immunogen	Human RBM14
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	RNA binding motif protein 14;RBM14;COAA;DKFZp779J0927;MGC15912;MGC31756;PSP2;SIP;SYTIP1; TMEM137 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RBM14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.