



RBM26 Antibody

Product Code	CSB-PA019417GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T8P6
Immunogen	Human RBM26
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	RNA binding motif protein 26;RBM26;ARRS2;C13orf10;FLJ20957;MGC133295;MGC133296;PRO1777;R P11-255E21.1;SE70-2;ZC3H17 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RBM26
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.