



REEP6 Antibody

Product Code	CSB-PA019545GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HR9
Immunogen	Human REEP6
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	receptor accessory protein 6;REEP6;C19orf32;DP1L1;FLJ25383;TB2L1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	REEP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.