



RIC8B Antibody

Product Code	CSB-PA019709GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVN3
Immunogen	Human RIC8B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	resistance to inhibitors of cholinesterase 8 homolog B (C. elegans);RIC8B;FLJ10620;MGC39476;RIC8;hSyn ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RIC8B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.