



RIOK1 Antibody

Product Code	CSB-PA019732GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRS2
Immunogen	Human RIOK1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	RIO kinase 1 (yeast);RIOK1;AD034;FLJ30006;MGC12903;bA288G3.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RIOK1
Target Details	This gene includes two alternatively spliced transcript variants, which encode different isoforms. The function of this gene has not been determined.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.