



RNF20 Antibody

Product Code	CSB-PA019869GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VTR2
Immunogen	Human RNF20
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ring finger protein 20;RNF20;BRE1;BRE1A;FLJ11189;FLJ20382;KIAA2779;MGC129667;MGC129668;hBRE1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RNF20
Target Details	This protein shares similarity with BRE1 of <i>S. cerevisiae</i> . Yeast BRE1 is a ubiquitin ligase required for the ubiquitination of histone H2B and the methylation of histone H3.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.