



RPGRIP1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA020100GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96KN7 |
| Immunogen | Human RPGRIP1 |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | retinitis pigmentosa GTPase regulator interacting protein 1; RPGRIP1; CORD9; DKFZp686P0897; LCA6; RG11; RGRIP; RPGRIP; RPGRIP1d ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | RPGRIP1 |
| Target Details | This gene encodes a photoreceptor protein that interacts with retinitis pigmentosa GTPase regulator protein and is a key component of cone and rod photoreceptor cells. Mutations in this gene lead to autosomal recessive congenital blindness. |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |