



# Recombinant Rabbit Voltage-dependent L-type calcium channel subunit beta-3 (CACNB3)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP004413RB  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P54286  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Oryctolagus cuniculus (Rabbit)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | <p>MYEDSYVPGF EDSEAGSADS YTSRPSLDSD VSLEEDRESA RREVESQAQQ<br/> QLERAKHKPV AFAVRTNSY CGVLDEECPV QGSGINFEAK DFLHIKEKYS<br/> NDWWIGRLVK EGGDIAFIPS PQRLESIRLK QEQRARRSGN PSSLSDIGNR<br/> RSPPPSLAKQ KQKQAEHVPP YDVVPSMRPV VLVGPSLKG Y EVTDMMQKAL<br/> FDFLKHFRFDG RISITRV TAD LSLAKRSVLN NPGKRTIIE R SSARSSIAEV<br/> QSEIERIFEL AKSLQLVLD ADTINHPAQL AKTSLAPIIV FVKVSSPKVL<br/> QRLIRSRGKS QMKHLTVQMM AYDKLVQCPP ESFDVILDEN QLEDACEHLA<br/> EYLEVYWRAT HHPAPGPGLL GPPSAIPGLQ SQQLLGERGE EHSPLERDSL<br/> MPSDEAWTGS SQRSSRHEE DYADAYQDLY QPHRQHTSGL<br/> PSANGHDPQD RLLAQDSEHN HNDNRWQRNR PWPKDSY</p> |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | CACNB3  |
| <b>Protein Names</b>     | Recommended name: Voltage-dependent L-type calcium channel subunit beta-3 Short name= CAB3 Alternative name(s): Calcium channel voltage-dependent subunit beta 3  |
| <b>Expression Region</b> | 1-477   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |