



Recombinant Methanococcus maripaludis Protein pelota homolog (pelA)

Product Code	CSB-BP017778MNR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6VI74
Product Type	Recombinant Protein
Immunogen Species	Methanococcus maripaludis (strain C7 / ATCC BAA-1331)
Purity	≥85% (SDS-PAGE)
Sequence	MKIIQEIPEK NIIKLIPENL DDLWHLNII QPYNAIYAVT ERRTEDKGDK LRADRGTKRK VFLGIKAEKI NFHEDFNRLR VSGKIIHAPE DIPIGSYHTI DIEPLLQVSV QKNWKKWDLA RLKDAEDSSK KPKVVVVIMD DSEADIFLVR EFGIKELASI KSGVSKKLDY KQNEQAKFSY YSDIINSISE FEGKILFAGP GFGKNNIQNY ISEKHKDLAP NVVVESANHT GKSGLSEILK SGIIDKIYGE ARISKETQII EKLLEEISKK GLAAYGIESV NNAMNYSIAID TLLLTDEYLR RNRRTIEELM NSVENINGSI LIISTEHDAG KQLKALGGIS ALLRFPIE
Source	Baculovirus
Target Names	pelA
Protein Names	Recommended name: Protein pelota homolog EC= 3.1.-.-
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.