



SAMD4A Antibody

Product Code	CSB-PA020684GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPU9
Immunogen	Human SAMD4A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	sterile alpha motif domain containing 4A;SAMD4A;DKFZp434H0350;KIAA1053;SAMD4;SMG;SMGA;Smaug;Smaug1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SAMD4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.