



SATB2 Antibody

Product Code	CSB-PA020720GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPW6
Immunogen	Human SATB2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;
Form	Liquid
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SATB homeobox 2;SATB2;FLJ21474;FLJ32076;KIAA1034;MGC119474;MGC119477 ;
Immunogen Species	Homo sapiens (Human)
Target Names	SATB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.