



SCARB1 Antibody

Product Code	CSB-PA020767GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WTV0
Immunogen	Human SCARB1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	scavenger receptor class B, member 1;SCARB1;CD36L1;CLA-1;CLA1;MGC138242;SR-BI;SRB1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SCARB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.