



SCLT1 Antibody

Product Code	CSB-PA020826GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NL6
Immunogen	Human SCLT1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	sodium channel and clathrin linker 1;SCLT1;CAP1A;FLJ30655;FLJ36042;MGC70542;hCAP-1A ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SCLT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.