



SCUBE3 Antibody

Product Code	CSB-PA020868GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IX30
Immunogen	Human SCUBE3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	signal peptide, CUB domain, EGF-like 3;SCUBE3;CEGF3;DKFZp686B09105;DKFZp686B1223;DKFZp686D20108;FLJ 34743 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SCUBE3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.