



SERF2 Antibody

Product Code	CSB-PA021042GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84101
Immunogen	Human SERF2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	small EDRK-rich factor 2;SERF2;4F5REL;FAM2C;FLJ20431;FLJ37527;FLJ38557;H4F5rel;HsT17089; MGC48826 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SERF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.