



SERPINA4 Antibody

Product Code	CSB-PA021060GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29622
Immunogen	Human SERPINA4
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 4;SERPINA4;KAL;KLST;KST;PI4;kallistatin ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SERPINA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.