



SERPINB2 Antibody

Product Code	CSB-PA021070GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05120
Immunogen	Human SERPINB2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	serpin peptidase inhibitor, clade B (ovalbumin), member 2; SERPINB2; HsT1201; PAI; PAI-2; PAI2; PLANH2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SERPINB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.