



SERPINE2 Antibody

Product Code	CSB-PA021082GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07093
Immunogen	Human SERPINE2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	serpin peptidase inhibitor, clade E (nexin, plasminogen activator inhibitor type 1), member 2;SERPINE2;GDN;PI7;PN1;PNI ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SERPINE2