



# Recombinant Human Segment polarity protein dishevelled homolog DVL-1-like (DVL1L1)

<b>Product Code</b>	CSB-YP007285HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P54792
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAETKIIYHM DEEETPYLVK LPVAPERVTL ADFKNVLSNR PVHAYKFFFK SMDQDFGVVK EEIFDDNAKL PCFNDRVSW LVLVEGAHSD AGSQGTDSHT DLPPPLERTG GIGDSRSPSF QPDVASSRDG MDNETGTESM VSHRRDRARR RNREEAARTN GHPRGDRRD VGLPPDSAST ALSSELESSS FVDSDEDDST SRLSSSTEQS TSSRLIRKHK RRRRKQRLRQ ADRASSFSSM TDSTMSLNII TVTLNMERHH FLGICIVGQS NDRGDGGIYI GSIMKGGAVA ADGRIEPGDM LLQVNDVNFENMSNDDAVRV LREIVSQTGP ISLTVAKCWD PTPRSYFTVP RPDPVRPIDP AAWLSHTAAL TGALPRPQLE EAPLTVESDM NTVVRVMQLP DSGLEIRDRM WLKITIANAV IGADVVDWLY THVEGFKERR EARKYASSLL KHGFLRHTVN KITFSEQCYY VFGDLCNLA TLNLNSGSSG TSDQDTLAPL PHPAAPWPLG QGYPYQYPGP PPCFPPAYQD PGFSYGS GST GSQQSEGSKS SGSTRNLRP PACEKERRAA GSGDSDSESD HTAPSGVGSS WRERPADQLS RGSSPRSQAS SYAPGLPPPH PTTKAYTVVG GPPGGPPVRE LAAVPELTG SRQSFQKAMG NPCEFFVDIM
<b>Source</b>	Yeast
<b>Target Names</b>	DVL1P1
<b>Protein Names</b>	Recommended name: Segment polarity protein dishevelled homolog DVL-1-like Short name= Dishevelled-1-like Alternative name(s): DSH homolog 1-like
<b>Expression Region</b>	1-670
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.