



SET Antibody

Product Code	CSB-PA021099GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01105
Immunogen	Human SET
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SET nuclear oncogene;SET;2PP2A;I2PP2A;IGAAD;IPP2A2;PHAPII;TAF-I;TAF-IBETA ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SET
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.