



SFPQ Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA021136GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P23246 |
| Immunogen | Human SFPQ |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | splicing factor proline/glutamine-rich (polypyrimidine tract binding protein associated);SFPQ;POMP100;PSF ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SFPQ |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |