



AKAP17A Antibody

Product Code	CSB-PA021150GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02040
Immunogen	Human SFRS17A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	splicing factor, arginine/serine-rich 17A; SFRS17A; 721P; CCDC133; CXYorf3; DXYS155E; MGC125365; MGC125366; MGC39904; XE7; XE7Y ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKAP17A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.