



SH2D4A Antibody

Product Code	CSB-PA021214GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H788
Immunogen	Human SH2D4A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SH2 domain containing 4A; SH2D4A; FLJ20967; SH2A ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SH2D4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.