



SH3D19 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA021228GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q5HYK7 |
| Immunogen | Human SH3D19 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IF |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | SH3 domain containing 19;SH3D19;EBP;EVE1;Kryn;MGC105136;MGC118910;MGC118911;MGC118912;MGC118913;SH3P19 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SH3D19 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |