



SH3GL3 Antibody

Product Code	CSB-PA021232GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99963
Immunogen	Human SH3GL3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SH3-domain GRB2-like 3;SH3GL3;CNSA3;EEN-2B-L3;EEN-B2;HsT19371;SH3D2C;SH3P13 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SH3GL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.