



SHC4 Antibody

Product Code	CSB-PA021257GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6S5L8
Immunogen	Human SHC4
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	SHC (Src homology 2 domain containing) family, member 4; SHC4; MGC34023; RaLP; SHCD ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SHC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.