



SLAMF7 Antibody

Product Code	CSB-PA021369GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ25
Immunogen	Human SLAMF7
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	SLAM family member 7; SLAMF7; 19A; CD319; CRACC; CS1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLAMF7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.