



SLC18A1 Antibody

Product Code	CSB-PA021426GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54219
Immunogen	Human SLC18A1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 18 (vesicular monoamine), member 1;SLC18A1;CGAT;VAT1;VMAT1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC18A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.