



SLC1A1 Antibody

Product Code	CSB-PA021432GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43005
Immunogen	Human SLC1A1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	solute carrier family 1 (neuronal/epithelial high affinity glutamate transporter, system Xag), member 1;SLC1A1;EAAC1;EAAT3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC1A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.