



SLC1A6 Antibody

Product Code	CSB-PA021437GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48664
Immunogen	Human SLC1A6
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	solute carrier family 1 (high affinity aspartate/glutamate transporter), member 6;SLC1A6;EAAT4;MGC33092;MGC43671 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC1A6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.