



SLC22A15 Antibody

Product Code	CSB-PA021448GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZD6
Immunogen	Human SLC22A15
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 22, member 15; SLC22A15; DKFZp761G0313; FLIPT1; PRO34686 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC22A15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.