



SLC25A1 Antibody

Product Code	CSB-PA021475GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53007
Immunogen	Human SLC25A1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 25 (mitochondrial carrier; citrate transporter), member 1;SLC25A1;CTP;SLC20A3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC25A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.