



SLC25A3 Antibody

Product Code	CSB-PA021499GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00325
Immunogen	Human SLC25A3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20?, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	solute carrier family 25 (mitochondrial carrier; phosphate carrier), member 3;SLC25A3;OK/SW-cl.48;PHC ;
Immunogen Species	Homo sapiens (Human)
Target Names	SLC25A3