



# SLC25A3 Antibody

<b>Product Code</b>	CSB-PA021499GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q00325
<b>Immunogen</b>	Human SLC25A3
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA,IHC
<b>Form</b>	Liquid
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	solute carrier family 25 (mitochondrial carrier; phosphate carrier), member 3;SLC25A3;OK/SW-cl.48;PHC ;
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SLC25A3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.