



Recombinant *Klebsiella pneumoniae* subsp. *pneumoniae* Lipoyl synthase (lipA)

Product Code	CSB-YP012927KAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6T689
Product Type	Recombinant Protein
Immunogen Species	<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> (strain ATCC 700721 / MGH 78578)
Purity	>85% (SDS-PAGE)
Sequence	MSKPIVMERG VKYRDADKMA LIPVKNVATE REALLRKPEW MKIKLPADSS RIQGIKAAMR KNGLHSVCEE ASCPNLAECF NHGTATFMIL GAICTRRCPF CDVAHGRPVA PDANEPQKLA QTIADMGLRY VVTVSVDRDD LRDGGAQHFA DCISAIREKN PSIKIETLVP DFRGRMDRAL DILTVTPPDV FNHNLENVPR LYRQVRPGAD YNWSLKLLER FKEAHPEIPT KSGLMVGLGE TNDEIIEVMR DLRRHGVTML TLGQYLQPSR HHLPVQRYVS PEEFEEMKAE AMAMGFTHAA CGPFVRSYH ADLQAKGMEV K
Source	Yeast
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.