



SLC27A4 Antibody

Product Code	CSB-PA021536GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1M0
Immunogen	Human SLC27A4
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	solute carrier family 27 (fatty acid transporter), member 4; SLC27A4; ACSVL4; FATP4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC27A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.