



# SLC30A7 Antibody

<b>Product Code</b>	CSB-PA021572GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8NEW0
<b>Immunogen</b>	Human SLC30A7
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC,IF
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	solute carrier family 30 (zinc transporter), member 7;SLC30A7;DKFZp686M0368;ZNT7;ZnT-7;ZnTL2 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SLC30A7
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.