



SLC37A4 Antibody

Product Code	CSB-PA021613GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43826
Immunogen	Human SLC37A4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 37 (glucose-6-phosphate transporter), member 4;SLC37A4;G6PT1;G6PT2;G6PT3;GSD1b;GSD1c;GSD1d;MGC15729;PRO0685;TRG19 ;; G6PT
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC37A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.