



SLC39A7 Antibody

Product Code	CSB-PA021636GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92504
Immunogen	Human SLC39A7
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 39 (zinc transporter), member 7;SLC39A7;D6S115E;D6S2244E;H2-KE4;HKE4;KE4;RING5;ZIP7 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC39A7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.