



SLC7A4 Antibody

Product Code	CSB-PA021716GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43246
Immunogen	Human SLC7A4
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 7 (cationic amino acid transporter, y+ system), member 4;SLC7A4;CAT-4;CAT4;HCAT3;MGC129976;MGC129977;VH ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC7A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.