



SMAD9 Antibody

Product Code	CSB-PA021794GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15198
Immunogen	Human SMAD9
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SMAD family member 9;SMAD9;MADH6;MADH9;SMAD8A;SMAD8B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SMAD9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.