



SPARC Antibody

Product Code	CSB-PA022486GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09486
Immunogen	Human SPARC
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Pig
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	secreted protein, acidic, cysteine-rich (osteonectin); SPARC; ON ;Osteonectin; BM-40
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SPARC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.