



# Recombinant *Sinorhizobium medicae* Argininosuccinate synthase (argG)

<b>Product Code</b>	CSB-EP002234STV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A6UE09
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Sinorhizobium medicae</i> (strain WSM419) ( <i>Ensifer medicae</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASHKDVKKV VLAYSGGLDT SIILKWLQTE LGAEVVTFTA DLGQGEELEP ARKKAEMLGI KEIYIEDVRE EFKVDFVFPF FRANAVYEGV YLLGTSIARP LISKHLIDIA RKTGADIAH GATGKGNDQV RFELSAYALN PDIKIIAPWR DWSFKSRTDL LEFAEKHQIP VAKDKKGEAP FSVDANLLHS SSEGKVLEDP AQEAPEYVHM RTISPEAAPD KATIIKVGFE RGDAVSIDGV RMSAATLLAK LNEYGRDNGI GRLDLVENRF VGMKSRGVYE TPGGTILLSA HRAIESITLD RGAHLKDEL MPRYAELIYY GFWFSPEREM LQAAIDKSQE HVEGEVTLKL YKGNMVMVGR ESGKSLYSDK LVTFEDDQGA YDQKDAAGFI KLNALRLRTL AARNR
<b>Source</b>	E.coli
<b>Target Names</b>	argG
<b>Protein Names</b>	Recommended name: Argininosuccinate synthase EC= 6.3.4.5 Alternative name(s): Citrulline--aspartate ligase
<b>Expression Region</b>	1-405
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.