



SPN Antibody

Product Code	CSB-PA022590GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16150
Immunogen	Human SPN
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	sialoporphin; SPN; CD43; GPL115; LSN ; GALGP; Leukosialin; Leukocyte sialoglycoprotein; Galactoglycoprotein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SPN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.