



SRFBP1 Antibody

Product Code	CSB-PA022660GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEF9
Immunogen	Human SRFBP1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	serum response factor binding protein 1;SRFBP1;BUD22;FLJ25286;P49;Rlb1;STRAP;p49/STRAP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SRFBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.