



# Recombinant *Saccharomyces cerevisiae* Eukaryotic translation initiation factor 3 subunit I (TIF34)

<b>Product Code</b>	CSB-EP007538STA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A6ZMK5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKAIKLTGHE RPLTQVKY NK EGDLLFSCSK DSSASVWYSL NGERLGTLDG HTGTIWSIDV DCFTKYCVTG SADYSIKLWD VSNGQCVATW KSPVPVKRVE FSPCGNYFLA ILDNVMKNPG SINIYEIERD SATHELTKVS EEEPIHKIITH EGLDAATVAG WSTKGKYIIA GHKDGKISKY DVSNNYEYVD SIDLHEKSIS DMQFSPDLTY FITSSRDTNS FLVDVSTLQV LKKYETDCPL NTAVITPLKE FIILGGGQEA KDVTTTSANE GKFEARFYHK IFEEEIGRVQ GHFGPLNTVA ISPQGTSYAS GGEDGFIRLH HFEKSYFDFK YDVEKAAEAK EHMQEAN
<b>Source</b>	E.coli
<b>Target Names</b>	TIF34
<b>Protein Names</b>	Recommended name: Eukaryotic translation initiation factor 3 subunit I Short name= eIF3i Alternative name(s): Eukaryotic translation initiation factor 3 39 kDa subunit homolog Short name= eIF-3 39 kDa subunit homolog
<b>Expression Region</b>	1-347
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.